

MODALITIES OVERVIEW

There are several factors to consider when choosing the right technology for your product. Understanding the variables of each technology and its effect on different materials should be considered as early in the product design process as possible.

RADIATION



ELECTRON BEAM

E-beam irradiation exposes products to high-energy electrons



TECHNOLOGY VARIABLES

- ORIENTATION TO BEAM
- BEAM ENERGY
- BEAM POWER
- PRODUCT CONVEYANCE
- PROCESSING TIME



PRODUCT CONSIDERATIONS

MATERIAL COMPATIBILITY
Qualified at maximal dose received during routine process

DENSITY
Moderate to low penetration on dense products



POSSIBLE EFFECTS

Reduced (potential) material impact effects than gamma and X-ray



EXPOSURE TIME TO RADIATION FIELD

SECONDS TO MINUTES
Time varies based on dose requirements

ISO 11137
ISO 13485



GAMMA IRRADIATION

Gamma irradiation exposes products to a Cobalt 60 radiation field



TECHNOLOGY VARIABLES

- PROCESSING TIME
- ISOTOPE ACTIVITY
- PRODUCT CONVEYANCE



PRODUCT CONSIDERATIONS

MATERIAL COMPATIBILITY
Qualified at maximal dose received during routine process

DENSITY
Good penetration on dense products



POSSIBLE EFFECTS

On certain materials, embrittlement, discoloration, change in viscosity due to irradiation



EXPOSURE TIME TO RADIATION FIELD

HOURS
Time varies based on dose requirements

ISO 11137
ISO 13485



X-RAY

X-ray irradiation uses ionizing energy from high-powered electron beam accelerators



TECHNOLOGY VARIABLES

- PRODUCT CONVEYANCE
- PROCESSING TIME



PRODUCT CONSIDERATIONS

MATERIAL COMPATIBILITY
Qualified at maximal dose received during routine process

DENSITY
Excellent penetration on dense products



POSSIBLE EFFECTS

Similar to gamma, but less pronounced due to improved dose uniformity and processing time



EXPOSURE TIME TO RADIATION FIELD

MINUTES
Time varies based on dose requirements

ISO 11137
ISO 13485

GAS



ETHYLENE OXIDE

Ethylene oxide sterilization uses a 3-stage gas process that includes pre-conditioning, sterilization, and post-conditioning (aeration)



TECHNOLOGY VARIABLES

- PROCESSING TIME
- TEMPERATURE
- HUMIDITY
- EO CONCENTRATION
- 100% EO OR POSITIVE PRESSURE PROCESS
- CHAMBER SIZE
- SINGLE CHAMBER OR MULTI-CHAMBER



PRODUCT CONSIDERATIONS

MATERIAL COMPATIBILITY
Very few compatibility concerns

Liquids generally not recommended

Packaging must be breathable

Sources of ignition (e.g. batteries, moving parts) require additional assessment

DENSITY
May impact distribution of some technology variables and therefore exposure time

GAS PATHWAY
Relies on vacuum to remove non-condensable gas and must reach surface requiring sterilization



POSSIBLE EFFECTS

Due to accumulation of EO residuals products may require prolonged aeration



EXPOSURE TIME

HOURS TO DAYS
Time varies based on validation approach and time required for degassing

ISO 11135
ISO 13485



VAPORIZED HYDROGEN PEROXIDE

VHP is a deep vacuum, low temperature 3-phase vapor process that includes conditioning, sterilization, and post-conditioning



TECHNOLOGY VARIABLES

- PROCESSING TIME
- TEMPERATURE
- HUMIDITY
- PRESSURE (VACUUM)
- VHP CONCENTRATION



PRODUCT CONSIDERATIONS

MATERIAL COMPATIBILITY
Compatible with most materials

Cellulose-based materials and highly absorbent materials are not compatible

Packaging must be breathable

VAPOR PATHWAY
Relies on vacuum to remove non-condensable vapor and must reach surface requiring sterilization

DENSITY
Surface sterilization process (limited penetration capabilities)



POSSIBLE EFFECTS

VHP breaks down safely into water and oxygen

Low residual levels

No known oxidation or discoloration effects



EXPOSURE TIME

HOURS
Time varies based on cycle requirements

ISO 22441
ISO 13485



MOIST HEAT (STEAM)

Steam autoclaving is a 3-stage high-temperature gas process that includes pre-conditioning, sterilization, and post-conditioning



TECHNOLOGY VARIABLES

- PROCESSING TIME
- TEMPERATURE
- SUBATMOSPHERIC PRESSURE
Pre-condition and post-condition
- SUPERATMOSPHERIC PRESSURE
Exposure phase



PRODUCT CONSIDERATIONS

LIQUID PRODUCTS
Can effectively sterilize contained liquids and gels (syringes, vials, bottles, etc.)

MATERIAL COMPATIBILITY
Chemically compatible with most materials

Ductility and melting of solids need to be considered due to high temperatures

GAS PATHWAY
Relies on vacuum to remove non-condensable gas (air) to achieve sterilant (steam) contact with internal surfaces



POSSIBLE EFFECTS

Damage to heat-sensitive, pressure-sensitive, or moisture-sensitive equipment

Melting, distortion, and potential corrosion

Water is the only residual / by-product



EXPOSURE TIME

MINUTES TO HOURS
Time varies based on product temperature limitations and validation approach

ISO 17665
ISO 13485